

Nurse Call Door Light



- **3 Color** can be used to distinguish patient from regular calls, advance calls, bathroom calls, and entering nursing status.

Working Voltage	DC18V
Installation Method	Wall Mounted
Body Size	Ø 102mm

The diagram illustrates the wiring for three rooms (Room 1, Room 2, Room 3) and a central server management area. Each room contains a bed extension, a door light, and an alarm button. The wiring connects these devices to a central POE 16 switch and a Server Management unit. Power adapters (18-24V) are used for the bed extensions and alarm buttons. Network cables connect the devices to the switch, and power cables connect them to the adapters. The diagram also shows a 'Door light' and 'Alarm Button' in Room 4, which are connected to the same network and power system.